



Falasept

Surface disinfectant for all common surfaces; liquid concentrate

Packaging	Art.-No.	per Box	per Palett
1 l PE-Bottle	3744	10	360
10 l PE-Canister	3745	-	60

Properties:

Suitable for the disinfection and cleaning of floorings and other waterproof surfaces. Falasept is virucidal, yeasticidal and to a limited extent virucidal. Effective against all enveloped viruses (corona, influenza, HBV etc.) and against rotavirus and norovirus.

Falasept is pH-neutral and has a high material compatibility so it is suitable for use on plexiglass, anodized aluminum, linoleum, PVC and other surfaces. Free from fragrance and colouring agents! Dosage according to application (see table).

Field of Application:

For application in commercial areas including food business, e. g. communal facilities, nursing homes, gyms, solariums, food processing facilities, canteen kitchens and others.

Application:

Remove heavy contaminations from the surface before use. Always make fresh diluted solution before use (see dosage table). Application by wiping with a mop or a cloth on surfaces. Let soak according to instruction (please note application time!). The surface needs to be wet for the entire application time. Let dry or rinse with clear water. Disinfected surfaces that come in contact with food products need to be rinsed with drinking water. The product must not be used in combination with other biocides or cleaning agents.

Working temperature: 5 – 25 °C


Ingredients:

5-15% cationic surfactants. Active substance in 100 g product: 3,3 g Didecylmethylammoniumchloride, 3,3 g Alkyl (C12-16) Dimethylbenzyl ammonium chloride, 3,3 g C12-14-Alkyl[(ethylphenyl)methyl]dimethyl-, chloride.

Productcode: GG 33

UFI: D5F0-U9WY-Y00J-900G

This product is intended for commercial use only. Safety data sheet available on request or at www.fala.de. For safety instructions pay attention to product label or safety data sheet.




Efficacy against	Concentration	Application time
Enveloped viruses, limited virucidity (RKI)	3,0 % 7,5 %	15 min. 1 min.
Non-enveloped viruses, virucidity against rota virus (RKI/DVV) noro virus (EN 14476)	0,25 % 6,0 %	1 min. 30 min.
Bactericidal, low load (EN 13697)	2,0 %	5 min.
Bactericidal, high load (EN 1276)	3,0 %	5 min.
Yeasticidal, effective against yeasts (EN 13697)	1,5 %	15 min.

Frequency of application: 1 – 2x per day or according to the current hygienic standards.

BAuA-Reg.-No.: N-94309, N-94310

Use biocides safely. Always read the label and product information before use.

The information contained herein corresponds at the time of publication to our best present-day knowledge. It is the responsibility of the user to ensure that the product is suitable for the intended application.



Disinfectant



Applicable for surface disinfection



Also suitable for food processing facilities